

Winemaking Analysis

Why and What to Test

John Rodrigues



Today's Presentation!

- **Who We are and What We Do**
- **Why Analyze**
- **When to Analysis**

Home Wine Lab

- We are a local laboratory specializing in routine wine analysis at reasonable prices.

We offer the following:

- Free Advice
- Next Day Results
- Courier Service
- Drop Off Locations:
(Old Possum Brewery/ Shone Farm / Sonoma Mission Gardens)
- Society Discount (10%)



www.homewinelab.com

(707) 978-4109

info@homewinelab.com

How we check your wine?

We use traditional analytical methods for checking your wine.



Aeration-Oxidation-Distillation
Apparatus



Automated Spectrometer



Cash Still

Instrument	Analysis
Digital Density Meter	Brix
Benchtop pH Meter	pH
Autotitrator	Titratable Acidity
Automated Spectrometer	Residual Sugar, L-Malic Acid, Yeast Assimilable Nitrogen
Electric Ebulliometer	Alcohol (Ethanol)
Aeration-Oxidation Apparatus	Free and Total Sulfur Dioxide
Cash Still	Volatile Acidity

Why is Analysis Important?

Best tools you have for winemaking are your senses.



Analysis is the best way to confirm your senses are correct!

Reduce Risk of Spoilage: Preventing Unwanted Characteristics



Hydrogen Sulfide (H_2S) is the Rotten Egg Odor



Over Sulfuring



Brettanomyces: Band-Aid, Barnyard, Sweaty Saddle




Oxidation



Volatile Acidity (VA): Vinegar or Nail Polish

Providing Numbers for Reproducibility

- It's not hard to make wine, but to make wine that you enjoy can be difficult.
- Having knowledge about the chemical makeup of your wine can help you achieve the style you want.




Wine Analysis Chart

- Major Winemaking Stages in Chronological Order
- Recommend Analysis for Each Stage
- Reason for Analysis
- Recommend Parameters for Each Type of Analysis

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titrateable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO₂ Based on pH Total SO₂ White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO ₂ TSO ₂	Informs SO ₂ adjustments	Volatile Acidity < 0.70 g/L
Pre-Bottling	pH TA FSO ₂ TSO ₂ RS VA Alc	Final state of wine and SO ₂ adjustments	Alcohol White 12 -14 % Red 13 - 15 % YAN 150 - 400 ppm

Grape Ripening



Development of sugars, acids, aromatics and berry color and weight during ripening

Visual Inspection:

Color Change (Veraison), Skin Soften, Seeds Turn Green to Brown

Chemical:

Brix = Increasing; TA = Decreasing;

pH = Increasing

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titrateable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO2 Based on pH Total SO2 White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO2 TSO2	Informs SO2 adjustments	Volatile Acidity < 0.70 g/L Alcohol White 12 - 14 % Red 13 - 15 %
Pre-Bottling	pH TA FSO2 TSO2 RS VA Alc	Final state of wine and SO2 adjustments	YAN 150 - 400 ppm

Crush



Visual:
Clean (No Bunch Rot, Skin Breakage, Raisining)

Chemical:
Homogenous Sample (Brix, pH, TA)
Nutrient Level (YAN = Amino Nitrogen & Ammonia)

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titratable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO2 Based on pH Total SO2 White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO2 TSO2	Informs SO2 adjustments	Volatile Acidity < 0.70 g/L
Pre-Bottling	pH TA FSO2 TSO2 RS VA Alc	Final state of wine and SO2 adjustments	Alcohol White 12 - 14 % Red 13 - 15 % YAN 150 - 400 ppm

YAN

°Brix of must or juice	Target YAN concentration (mg/L)
21	200
23	250
25	300
27	350


Lack of Nitrogen:

Sluggish Fermentation, Off-Odors (H_2S), Stuck Fermentation

Too Much Nitrogen:

Production of Acetate Ester (Nail Polish Aroma), Leftover Food for Bacteria

Fermentation




Visual Inspection:

CO₂ Production (Rise in Cap) (Bubbles)

Chemical:

Brix = Decreasing; Temperature = Increasing (not above 90 F)



* Add Nutrients Between 18 to 13° Brix

Confirm Completed Fermentation



Visual Inspection:

Cap Dropping (Sinking Skins)

Chemical:

Brix Stop Moving (Same for 2-3 days in a row)

Temperature = Decreases

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titrateable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO ₂ Based on pH Total SO ₂ White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO ₂ TSO ₂	Informs SO ₂ adjustments	Volatile Acidity < 0.70 g/L
Pre-Bottling	pH TA FSO ₂ TSO ₂ RS VA Alc	Final state of wine and SO ₂ adjustments	Alcohol White 12 - 14 % Red 13 - 15 % YAN 150 - 400 ppm

Primary & Secondary Fermentations

**Primary Fermentation =
Glucose/ Fructose**


**Secondary Fermentation =
Malic Acid**

Ensure both fermentations reach
safe concentration levels to
prevent future problems.

Chemical:

Residual Sugar (G+F) = at or below 1 g/L

L-Malic Acid = at or below 0.30 g/L



Post-Fermentation



Visual Inspection:

Taste Like Wine (Not Sweet, Too Tart, Vinegar Aroma)

Chemical:

pH and TA = Stability/ Preference; VA = Spoilage;

Alcohol = Matches Intention (Benchmark)

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titrateable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO2 Based on pH Total SO2 White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO2 TSO2	Informs SO2 adjustments	Volatile Acidity < 0.70 g/L Alcohol White 12 -14 % Red 13 - 15 %
Pre-Bottling	pH TA FSO2 TSO2 RS VA Alc	Final state of wine and SO2 adjustments	YAN 150 - 400 ppm

Maturing: Barrel/ Carboy



Prevent spoilage by monitoring and maintaining effective SO₂ levels

Chemical:

Free & Total SO₂ = Maintaining Effective Levels

pH = Determines Target Free SO₂ Level

* Volatile Acidity = Monitor Occasionally

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25
Crush	Brix pH TA YAN	Shows needed adjustments	pH White 3.20 - 3.40 Red 3.40 - 3.70
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Titratable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Free SO ₂ Based on pH
Maturing: Barrel/Carboy	pH FSO ₂ TSO ₂	Informs SO ₂ adjustments	Total SO ₂ White < 100 ppm Red < 120 ppm
Pre-Bottling	pH TA FSO ₂ TSO ₂ RS VA Alc	Final state of wine and SO ₂ adjustments	Residual Sugar ≤ 1 g/L
			L-Malic Acid ≤ 0.30 g/L
			Volatile Acidity < 0.70 g/L
			Alcohol White 12 -14 % Red 13 - 15 %
			YAN 150 - 400 ppm


How Does SO₂ Work?

Free SO₂:

- Active form that contributes to wine's stability

Bound SO₂:

- Inactive form that doesn't have any antioxidant or antimicrobial activity



Total SO₂:

- Free and Bound SO₂ added together

SO₂ Chart

Free SO₂ = Active Part of SO₂

pH = Affects equilibrium of cells; lower pH means higher percentage of molecular SO₂ effectiveness

Molecular SO₂ = Only part that is antimicrobial; small part of Free SO₂

***Sensory Effects: Too Much!!!**

- Chemical Taste (> 100 ppm Total)
- Burnt Matchstick Smell
 - > 0.8 to 1 ppm molecular SO₂

pH	White Wine 0.8 ppm molecular SO ₂	Red Wine 0.5 ppm molecular SO ₂
3	13	8
3.05	15	9
3.1	16	10
3.15	19	12
3.2	21	13
3.25	23	15
3.3	26	16
3.35	29	18
3.4	32	20
3.45	37	23
3.5	40	25
3.55	46	29
3.6	50	31
3.65	57	36
3.7	63	39
3.75	72	45
3.8	79	49
3.85	91	57
3.9	99	62

▲ FIGURE a Amount of “free SO₂” Required to Obtain the Correct Amount of Molecular SO₂ Concentration as a Function of pH in White and Red Wine.

Pre-Bottling



Last chance to change anything!

Reproducibility: reference point for future vintages

Final SO₂ Adjustment: Add 5 ppm Extra

Winemaking Stage	Analysis	Reason	Analysis Parameters
Grape Ripening	Brix pH TA	Helps determine ripeness	Brix White 21 - 24 Red 23 - 25 pH White 3.20 - 3.40 Red 3.40 - 3.70
Crush	Brix pH TA YAN	Shows needed adjustments	Titratable Acidity White 6 - 9 g/L Red 6 - 7 g/L
Confirm Completed Fermentation	RS ML	Indicates completed fermentation	Free SO₂ Based on pH Total SO₂ White < 100 ppm Red < 120 ppm
Post-Fermentation	pH TA Alc VA	Adjustments and potential spoilage	Residual Sugar ≤ 1 g/L L-Malic Acid ≤ 0.30 g/L
Maturing: Barrel/Carboy	pH FSO ₂ TSO ₂	Informs SO ₂ adjustments	Volatile Acidity < 0.70 g/L Alcohol White 12 -14 % Red 13 - 15 %
Pre-Bottling	pH TA FSO ₂ TSO ₂ RS VA Alc	Final state of wine and SO ₂ adjustments	YAN 150 - 400 ppm

Thank You!

10% Discount on All Analysis for Winemaking Societies